

GLOVE SELECTION CHART

The following guide is a general guide for glove selection in relation to chemicals handled. The information presented here is believed to be accurate; however, we cannot guarantee its accuracy. Many factors affect the breakthrough times of glove materials including, but not limited to:

1. Thickness of glove material
2. Concentration of the chemical worked with
3. Amount of chemical the glove comes in contact with
4. Length of time which the glove is exposed to the chemical
5. Temperature at which the work is done
6. Possibility of abrasion or puncture.

Some Common Sense Rules for Glove Use

- Select gloves which are resistant to the chemicals you may be exposed to. Consult the relevant Safety Data Sheet (SDS) which may recommend a particular glove material.
- Select gloves of the correct size and fitting; gloves that are too small are uncomfortable and may tear whereas overlarge gloves may interfere with dexterity. In some cases, such as use of HF, it may be advisable to select gloves that can be removed very rapidly in an emergency.
- Before use, check gloves (even new ones) for physical damage such as tears or pin holes and for previous chemical damage: this is especially important when dealing with dangerous materials such as HF.
- When working, it may be advisable to wash the external surface of the gloves frequently with water.
- Some gloves, especially lightweight disposables, may be flammable: keep hands well away from naked flames or other high temperature heat sources.
- When removing gloves, do so in a way that avoids the contaminated exterior contacting the skin.
- Wash hands after removing gloves.
- Dispose of contaminated gloves properly.
- Do not attempt to re-use disposable gloves.
- Never wear possibly contaminated gloves outside of the laboratory or to handle telephones, computer keyboards, *etc.*

Selection Key:

- 4** Excellent, breakthrough times generally greater than 8 hours.
- 3** Good, breakthrough times generally greater than 4 hours.
- 2** Fair, breakthrough times generally greater than 1 hour.

- 1** Not Recommended, breakthrough times generally less than 1 hour.
- ?** Not Tested or Information unknown. Use known tested glove type.

GLOVE SELECTION GUIDE

	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Chemical						
Organic Acids						
Acetic acid	2	3	4	2	1	4
Formic acid	2	3	4	3	2	2
Lactic Acid	4	4	4	3	4	4
Maleic acid	3	3	2	3	3	4
Oxalic acid	4	4	4	4	4	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Inorganic acids						
Chromic acid up to 70%	1	1	4	3	3	4
Hydrochloric acid up to 37%	3	3	4	3	3	3
Hydrofluoric acid up to 70%	2	2	3	1	1	?
Nitric acid 70+ %	?	1	2	?	1	4
Perchloric acid up to 70%	4	4	3	4	4	4
Phosphoric acid 70+ %	4	4	4	4	4	4
Sulfuric acid 70+ %	1	2	4	2	1	2

	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Alkalis						
Ammonium hydroxide up to 70%	1	3	4	2	3	?
Potassium hydroxide up to 70%	4	4	4	4	4	4
Sodium hydroxide 70+ %	4	4	4	4	3	3
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Salt Solutions						
Ammonium nitrate	4	4	4	4	4	4
Calcium hypochlorite	1	3	4	4	3	4
Ferric chloride	4	4	4	4	4	4
Mercuric chloride	3	3	4	3	3	4
Potassium cyanide	4	4	4	4	4	4
Potassium dichromate	4	4	4	4	4	4
Potassium permanganate	4	4	?	4	4	?
Sodium cyanide	4	4	4	4	4	4
Sodium thiosulfate	4	4	4	4	4	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Aromatic hydrocarbons						
Benzene	1	1	1	1	1	3
Gasoline	1	1	1	1	4	4
Naphthalene	1	1	1	1	4	4

Toluene	1	1	1	1	1	4
Xylene	1	1	1	1	1	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Aliphatic hydrocarbons						
Diesel fuel	1	2	1	2	3	4
Hexanes	1	1	1	1	4	4
Kerosene	1	3	1	3	4	4
Naphtha	1	2	1	3	4	4
Pentane	1	1	1	1	3	4
Petroleum ether	1	1	1	2	3	4
Turpentine	1	1	1	1	2	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Halogenated hydrocarbons						
Carbon tetrachloride	1	1	1	1	1	4
Chloroform	1	1	1	1	1	4
Methylene chloride	1	1	1	1	2	3
Polychlorinated biphenyls (PCB's)	1	4	4	?	2	4
Perchloroethylene	1	1	1	1	2	4
Trichloroethylene	1	1	1	1	1	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □

Esters						
Ethyl acetate	1	1	3	1	1	1
Butyl acetate	1	1	2	1	1	1
Methyl acetate	1	1	4	1	1	1
Isobutyl acrylate	1	1	4	1	1	1
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Ethers/Glycols						
Diethyl ether	1	2	1	1	2	1
Ethylene glycol	1	2	4	1	2	4
Isopropyl ether	1	2	1	1	3	1
Propylene glycol	?	3	3	2	2	?
Tetrahydrofuran	1	1	2	1	1	1
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □
Aldehydes						
Acetaldehyde	1	1	4	1	1	1
Acrolein	1	1	4	1	1	1
Benzaldehyde	1	1	4	1	1	3
Butyraldehyde	1	1	4	1	1	1
Formaldehyde	1	2	4	2	4	4
Glutaraldehyde	?	4	4	2	?	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton □

Ketones						
Acetone	1	1	4	1	1	1
Diisobutyl ketone	1	1	2	1	1	2
Methyl ethyl ketone	1	1	4	1	1	1
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton□
Alcohols						
Allyl alcohol	1	1	4	1	4	3
Butyl alcohol	1	3	4	2	3	4
Ethyl alcohol	1	2	4	1	3	4
Isopropyl alcohol	1	3	4	2	4	4
Methyl alcohol	1	1	4	1	1	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton□
Amines						
Aniline	1	1	4	1	1	2
Ethanolamine	2	4	4	3	4	4
Ethylamine	1	2	4	1	1	1
Methylamine	1	3	4	2	4	4
Triethanolamine	1	1	4	1	4	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton□
Elements						
Bromine	1	2	1	?	1	4

Chlorine aqueous	?	1	2	?	1	4
Iodine	?	1	3	?	3	4
Mercury	?	4	4	?	4	4
	Natural Rubber	Neoprene	Butyl	PVC	Nitrile	Viton [®]
Miscellaneous						
Acetic anhydride	1	2	4	1	1	1
Acetonitrile	1	1	4	1	1	1
Acrylamide	1	1	3	1	2	3
Carbon disulfide	1	1	1	1	1	4
Cresols	1	3	4	?	2	4
Cutting fluid	?	2	?	2	3	?
Dimethyl sulfoxide	1	4	4	1	1	1
Hydraulic oil	?	?	1	2	3	?
Hydrazine	2	4	4	4	4	1
Hydrogen Peroxide	4	2	4	3	4	4
Lubricating oil	3	3	?	?	4	3
Malathion	?	3	1	?	3	?
Nitrobenzene	1	1	4	1	1	4
Phenol	1	3	2	1	1	4
Photo solutions	3	4	?	3	4	?
Picric acid	1	2	3	1	2	4
Pyridine	1	1	4	1	1	1

Viton is a registered trademark of DuPont Dow Elastomers.



University of Wyoming Safety
Web: www.uwyo.edu/safety/
Phone: (307) 766-3277
Email: uwehs@uwyo.edu
Regulated Materials Management Center
Phone: (307)766-3698 Fax: (307)766-3699
Email: HAZMAT@uwyo.edu

This information was obtained from <http://www.allsafetyproducts.biz/page/74172>